# ANALYSE LIPID NANOPARTICLES AND OTHER DRUG DELIVERY VEHICLES ON A PARTICLE-BY-PARTICLE LEVEL



## THE EXOID: A POWERFUL TOOL FOR PARTICLE ANALYSIS IN NANOMEDICINE

- Simultaneously measure size and concentration, or size and zeta potential of lipid nanoparticles and other drug delivery vehicles in solution
- The single-particle nature of Tunable Resistive Pulse Sensing (TRPS) enables unmatched resolution, highly suited to analysing polydisperse samples
- Accurate particle size distribution analysis lends the Exoid to stability studies



#### PRECISE PARTICLE SIZE DISTRIBUTION ANALYSIS

Figure 1. Size of empty LNPs and mRNA-loaded LNPs. (A) Representative size distribution graph generated using the Exoid. (B) Median size ± interquartile range of empty LNPs (n=4) and mRNA-loaded LNPs (n=4), compared using the Mann Whitney U test, \*p<0.05





#### COMPARE SAMPLES WITH CONFIDENCE

**Figure 2**. Size of fresh and freeze-thawed mRNA-loaded LNPs. (A) Representative size distribution graph generated using the Exoid. (B) Median size and interquartile range of fresh mRNA-loaded LNPs (n=4) and freeze-thawed mRNA-loaded LNPs (n=3), compared using the Mann-Whitney U test, \*p<0.05.



### MEASURE SIZE AND ZETA POTENTIAL



Figure 3. Zeta potential of negatively-charged empty and mRNA-loaded lipid nanoparticles

4

## THE EXOID BRINGS SINGLE-PARTICLE CHARACTERISATION TO NANOMEDICINE

Assess polydispersity

Capture the heterogeneity of particle size as part of your quality control process.

Optimise process development

Identify subtle changes to size, concentration or zeta potential that might reflect aggregation during development or storage.

▶ Get deeper insights

Obtain high-resolution, single-particle insights rather than bulk estimates.

0

5



200

## A STANDARDISED APPROACH TO NANOPARTICLE ANALYSIS



Reproducibility enabled through the use of standardised NISTtraceable calibration particles

- Clean user interface
- Automated data processing

No reliance on prior knowledge of optical properties of particles or dispersant





# UNMATCHED ABILITY TO RESOLVE POLYDISPERSE SAMPLES



**Figure 4**. Comparison of Tunable Resistive Pulse Sensing, Nanoparticle Tracking Analysis, and Multi-Angle Dynamic Light Scattering (MADLS) measurements of a quadrimodal sample of polystyrene bead standards.



UNMATCHED RESOLUTION. SINGLE-PARTICLE INSIGHTS.

WWW.IZON.COM/TRPS/EXOID