
# © © MAGNETIC CONCENTRATION KIT

FAST, AUTOMATABLE EV CONCENTRATION






### A new age of EV concentration is here: the qEV Magnetic Concentration Kit

- Concentrate EVs that have been purified using qEV columns.
- Using functionalised magnetic particles, you can concentrate EVs more quickly and cleanly than before without the need for centrifugation or time-consuming protocols.

### A quick and automatable approach to concentrating EVs

Why spend 75 – 90 minutes using **other EV concentration methods** such as centrifugal concentration filters?

With the **qEV Magnetic Concentration Kit**, you can concentrate your EVs in just 20 minutes.



#### Concentrate efficiently with the qEV Magnetic Concentration Kit

- ▶ All-in-one system for concentrating intact EVs
- Uses Magnetic Nanotrap® Extracellular Vesicle Particles functionalised with chemical affinity baits
- Does not require a centrifuge, making the qEV Magnetic Concentration Kit amenable to automation

#### What is the qEV Magnetic Concentration Kit suitable for?

- Omics studies (e.g., proteomics, RNA sequencing)
- Western blot
- Functional studies
- ELISA



## CAPTURE AND CONCENTRATE YOUR ISOLATED EVS.

#### LEARN MORE NOW AT



