

INSTANT DOCUMENT ACCESS:



support.izon.com/getting-started-with-the-exoid

SCAN ME

Use the **QR code** or **URL** to access the latest user manuals and safety data sheets for the Exoid.



IMPORTANT PRE-USE INFORMATION

Installation and training

- ▶ If using the TRPS Training Kit, please ensure all components are stored correctly prior to use.
- ▶ For IQ/OQ installation: your support representative will be in touch to arrange installation and training.
- ▶ For non-IQ/OQ installation: follow the instructions detailed in the Exoid training manual, which can be followed alongside the Izon Academy training (academy.izon.com) for a comprehensive tutorial.

Support

- ▶ For additional assistance, visit support.izon.com
- ▶ Any further enquiries or repair/maintenance requests should be made to our support team via the online portal.
- ▶ When reporting Exoid issues, please ensure you provide the serial number of the Exoid.


Please do not dispose of Exoid packaging

If you need to return the Exoid and/or components, you must use the original packaging to prevent damage during shipment. Please retain all packaging for this purpose.



WARNING

Refer to the Exoid user manual for all warnings related to its use.



1. Exoid Instrument
2. Shielding Cap
3. Upper Fluid Cell
4. 10 M Ω Reference Cell

5. Lower Fluid Cell
6. 24 V Power Supply
7. Power Lead
8. USB-A to USB-C 3.0

