

## THE EXOID



Before unpacking and installing your Exoid there are are few important things to note:

- If you require an IQ/OQ (also known as 3Q) installation, please remove the Solution S and TKP200 particles from the Training Kit and store them at 4-8°C. Close the box and set it aside for when your certified Izon Science IQ/OQ executor is present.
- If you **do not** require an IQ/OQ installation you can unpack and set up the rest of your system according to the instructions in the Exoid User Manual which can be downloaded from support.izon.com/exoid-user-manual
- The User Manual has detailed instrument descriptions and installation information which you will need to read through in its entirety before operating the instrument. Instructions on downloading the appropriate software suites can be found in the Assembly and Setup Instructions section.
- For non-IQ/OQ installations, download the training manual from support.izon.com/exoid-training-manual
- For extra support, proceed to support.izon.com

www.izon.com Rev E