

GETTING STARTED



WELCOME TO THE IZON FAMILY

Before you can get started on unpacking and installing your exciting new system, there are a few important things to cover:

- ▶ If you require an IQ/OQ (also known as 3Q) installation, please remove the Solution S and TKP200 particles from the Training Kit and store them at 4-8°C. Close the box and set it aside for when your certified Izon Science IQ/OQ executor is present.
- ▶ If you **do not** require an IQ/OQ installation you can unpack and set up the rest of your system according to the instructions in the qNano User Manual (for non-IQ/OQ installations) which can be downloaded from support.izon.com/reference-material-qnano-user-manual
- ▶ For non-IQ/OQ installations, download the Izon Control Suite from support.izon.com/how-can-i-get-the-latest-software-release
- ▶ For non-IQ/OQ installations, download the training manual from support.izon.com/overview-qnano-training-programme